

BA⁻

Minerals	Stock Solutions (g/L)	Quantity (mL Stock/L Media)	Final Concentration (mM)
NaCl	250	33.3	142
MgCl ₂ · 6H ₂ O	200	3.3	3.3
KCl	50	3.3	2.14
MgSO ₄ · 7H ₂ O	350	3.3	4.62
CaCl ₂ · 2H ₂ O	50	3.3	0.99
Na ₃ -citrate	0.6	8.3	0.027
Na ₂ -EDTA · 2H ₂ O	0.1	8.3	0.002
Trace metal mix (A5 + Co)	See recipe below	1	-

Adjust to 900mL with mQ water and autoclave.

After cooling, add the following filter sterilized (0.2 µm) components to complete the medium:

Minerals	Stock Solutions (g/L)	Quantity (mL Stock/L Media)	Final Concentration (mM)
K ₂ HPO ₄ · 3H ₂ O	4	8.3	0.15
Na ₂ CO ₃	20	2.3	0.44
Fe-NH ₄ -citrate	6	0.8	-
Vitamin B12 (Cyanocobalamin)	0.02	1	-

For solid medium use 7g/L of agarose. Sterilize the agarose separately in 550 ml of milliQ water. In this case the mineral solution is filled up to 400 ml.

Trace metal mix A5 + Co:

Trace metals	Quantity g/L	Concentration in the final media (mM)
H ₃ BO ₃	2.86	0.047
MnCl ₂ · 4H ₂ O	1.81	0.009
ZnSO ₄ · 7H ₂ O	0.22	0.0007
Na ₂ MoO ₄ · 2H ₂ O	0.39	0.0016
CuSO ₄ · 5H ₂ O	0.08	0.0003
Co(NO ₃) ₂ · 6H ₂ O	0.05	0.0002